



MEDICLIC HE

RANGE OF TUNGSTEN SYRINGE SHIELDS FOR RADIOPHARMACEUTICAL SUBSTANCES

Compatibility: with all syringes on the market.

Simplicity: can be completely disassembled by the operator for decontamination or to change the glass.

Ergonomics: lead glass window with zoom effect and white lining for optimum viewing.

Rapidity: quick and easy syringe insertion/withdrawal system (1 click).



Characteristics

General	1 mL short	1 mL long	2 mL	2 mL "BD Emerald"	2.5 mL	3 mL short	3 mL long	5 mL short	5 mL long	10 mL short	10 mL long
Size planned for the syringe	70 mm	74 mm	52 mm	49 mm	58 mm	62 mm	68 mm	59 mm	62 mm	73 mm	82 mm
Radiation protection:	5 mm tungsten / 9 mm lead glass										
Package dimensions:	L 240 x D 180 x H 50 mm										
Package weight:	0.130 kg + weight of syringe shield										
Ref.:	11653	11651	11663	11666	11661	11693	11691	11673	11671	11683	11681

Consumables

Mediclic tungsten glass TEP 1 mL short and long

Ref. 11652

Mediclic tungsten glass TEP 2 mL and 2 mL "BD Emerald"

Ref. 11664

Mediclic tungsten glass TEP 2.5 mL

Ref. 11662

Mediclic tungsten glass TEP 3 mL short

Ref. 11694

Mediclic tungsten glass TEP 3 mL long

Ref. 11692

Mediclic tungsten glass TEP 5 mL short and long

Ref. 11672

Mediclic tungsten glass TEP 10 mL short

Ref. 11684

Mediclic tungsten glass TC 10 mL long

Ref. 11682



Download the product data sheet
medisystem.fr