

MacroMedics®
UniqFix™

The UniqFix is a multifunctional carbon fiber baseplate, offering countless possibilities for the immobilisation of extremities. Dedicated thermoplastic mask shapes have been designed for immobilisation of arms, hands, legs and feet in straight or bent position. The UniqFix can be indexed on the couchtop in two perpendicular directions, allowing for arm positioning extended from the couchtop.

An MR safe version of the UniqFix is available.

135100	UNIQCF	UniqFix positioning system for extremities, carbon fiber
135150	UNIQFMR	UniqFix positioning system for extremities, MR safe

