

Conventional Simulation Products:



X-SPOTS®

order code: 101

contents per box: 150

1.5mm pellet with easy pik-up tab

perfect for marking outer canthus and nodes



Y-SPOTS®

order code: 102

contents per box: 132

2.0mm pellet with easy pik-up tab

perfect for marking outer canthus and nodes



V-SPOTS®

order code: 103

contents per box: 130

2.5mm pellet with easy pik-up tab

perfect for marking port films




Z-SPOTS®

order code: 114

contents per box: 40

4.3mm pellet with easy pik-up tab

perfect for marking port films



radiopaque arrow



D-SPOTS® arrow

order code: 104


contents per box: 99

radiopaque arrow with easy pik-up tab

perfect for marking nodes and outer canthus



.8mm diameter radiopaque line



T-SPOTS® line

order code: 111

contents per box: 275cm

.8mm diameter radiopaque line

perfect for marking scars, nodes, field borders and match lines