

Pediatric

AttenuRad CT Breast Shields



Studies Show AttenuRad CT Shields Provide Significant IN-PLANE Radiation Protection

FL® F&L Medical Products
PH: (724) 845-7028
FAX: (724) 845-5439
1129 Industrial Park Road
Vandergrift, PA 15690

www.fandlmedicalproducts.com

Composition:

F&L AttenuRad CT Shields are:

- Non-Latex
- Lead-Free




Sizes:

| | | |
|-----------|---------|---------|
| * Neonate | AR14X5 | 10/pack |
| * 1 year | AR24X6 | 8/pack |
| * 5 year | AR28X8 | 6/pack |
| * 10 year | AR33X10 | 4/pack |

*All Pediatric Shields and Adult Breast Shields come with .635cm foam offset. Additional offset foam available.

Protection:

.035mm Pb equivalent



Courtney Coursey, MD, Donald Frush, MD, et al, Society of Pediatric Radiology, 2006 Convention, Montreal, Quebec, Canada - May 2006.

