

SAFE SIMPLE AND SECURE...



Lead Glass for Vial Protection

Developed for low-energy isotopes and beta protection

- **Lead glass offers clear 360° visibility**
- **Accommodates most vials**
- **Small mechanical dimensions, smallest on the market**
- **Low weight**

The high density (5.2g/cm3) Lead Glass Vial Shield protects the user from hand and eye exposure, when inspecting vials before aspiration, and the new vial glass offers complete 360° visibility at all times. The vial glass is shielded for low-energy isotopes. Tungsten parts in top and bottom ensure the small physical dimension and weight, which offers maximum ergonomical comfort.

The lead glass has a unique magnetic lift-off top for needle injection. The enclosed vial adaptors and memory foam ensure that the vial is always correctly positioned and aligned.

Low-Energy Lead Glass

The low-energy lead glass is offered in two different lead glass thicknesses, 3 mm and 6 mm Pb equivalent.

Beta Lead Glass

The lead glass for Zevalin offers maximum protection with lead glass protection of 22.5 mm and an internal diameter of 37 mm for acrylic insert.



Part no	Description
10-001	Low activity vial glass
10-002	Med/High activity lead glass
10-003	Beta shield
10-004	Tray for 2 vial glasses
10-005	Tray for 3 vial glasses
10-006	Extra vial adaptor
10-007	Memory foam disc set

Optional:

- The lead glass vial is prepared for a counterbalancer for maximum ergonomic use.
- Tray for 2 and 3 vial glasses.
- Electrical rocker for mixing vial kit, 1 or 2 arms.
- Optional extra plastic vial adaptors, special adaptors and colours on request.
- Memory foam disc (set of 20).
- In case of defects, the vial glass can be repaired.

Specification

Lead equivalency 5.2g/cm3

Accommodates vials from 5 thru 30ml

Low Activity Vial Shield

- Lead glass (thickness 12.5 mm) equivalent to 3 mm Pb
- Length 9.00 cm
- Diameter 6.50 cm
- Weight 1195g

Medium/High Activity

- Lead glass (thickness 20 mm) equivalent to 6 mm Pb
- Length 9.00 cm
- Diameter 7.50 cm
- Weight 1700g

Beta Shield

- Length 11.00 cm
- Diameter 9.30 cm
- Weight 2100g



Lead Glass Vial Shield