

Stretcher Trendelenburg hydraulic variable height



Ref. 100 1002

TECHNICAL SPECIFICATIONS

- Fully tested and compatible 3 Tesla MRI
- Height adjustable with pedals from 760 to 1000 mm for ergonomic working
- Inclination of 12% - 12%
- Maximum load: 180 kg
- Protective barriers retractable to facilitate patient access
- Infusion stand 2 branches
- 4 wheels Ø 150 mm antistatic (including 3 w/brakes and 1 directional)
- Forend lifting of 0 to 75 °
- Mattress thickness 50 mm, Fire-resistant coating M1
- Coating facilitates cleaning and disinfection

