

MACROMEDICS®

Breastboard LX™

The MACROMEDICS® Breastboard LX is a high-end breastboard, designed to meet the latest requirements in breast treatment and manufactured with use of state of the art 3-D carbon fiber production techniques. The ergonomic design is a result of close cooperation with- and competent input from the clinical staff with maximum focus on treatment set-up, dose delivery, efficiency and patient comfort. The Breastboard LX provides a coreless single layer carbon fiber treatment zone without grid to obtain low-attenuation, maximal homogeneity, rigidity and stability. The Breastboard LX is a modern, very light and ergonomic breast treatment device securing reproducible day to day patient set-up.

The Breastboard LX features:

- Rigid, homogeneous and low-attenuation IGRT treatment zone
- Stable and easy to operate arm supports providing 85° - 150° angles to the midplane and indexed width-adjustment for the upper arm + fully adjustable supports for the lower arm
- No inserts or obstructing arm supports in the treatment area
- Rigid and easy to handle elevation brackets providing breast board angles from 50° - 200° in 2.50 steps, the brackets fold in the bottom plate for easy handling
- Ergonomic and stable lower body support with precision movement, adjustable in length in 10 mm steps
- Compatible with MACROCAST™ 2-point head-, neck- and breast support thermoplastics
- Comes with a comfortable wide head support, the position of the head support can be adjusted
- The Breastboard LX is indexable and compatible with all actual couch tops for radiotherapy and X-ray settings.

Contact MACROMEDICS or your local distributor for detailed information or an on site product demo.

122500	MBLXI	Breastboard LX
--------	-------	----------------



macromedics

SOLUTIONS IN RADIOTHERAPY

MACROMEDICS BV
Kouwe Hoek 18
2741 PX Waddinxveen
The Netherlands

Tel: +31(0)182 389777
Fax: +31(0)182 389778
info@macromedics.com
www.macromedics.com