

MACROMEDICS®

Exafix IMRT head- and neck Baseplates™

Intensity Modulated Radiation Therapy requires high precision and reliable and accurate patient-set-up. MACROMEDICS® offers a IMRT baseplate specifically designed to fit these requirements.

The Exafix IMRT Baseplates™ are developed and produced with state of the art special 3-D molding- and production techniques resulting in user-friendly shapes with focus on the significant and important areas with respect for low tolerances in attenuation and dosimetry. The extended baseplate with ergonomic design with double curved surfaces supports the full upper body, offers a smooth transition in the lower back and guarantees a comfortable and unstressed patient set-up.

The cantilevered design offers a '360 degrees angle of freedom' around the region of interest.

The Carbon fiber Exafix IMRT Baseplates are available for several types of treatment couches, CT and PET-CT.

The MACROMEDICS Exafix IMRT Baseplate is a state of the art, unique, ultra light treatment device securing reproducible day-to-day patient set-up.

Exafix IMRT Baseplates for head and neck:

- Better target coverage
- Optimal positioning and approach
- Easily hooked onto cranial end of treatment couch
- Low tolerances in attenuation and dosimetry
- Uniformity on all treatment modalities

110110	EXF-IMRT/VAIGRT	Exafix IMRT Baseplate, fit for Varian IGRT couchtop
110130	EXF-IMRT/VAEX	Exafix-IMRT baseplate, fit for Varian Exact couch
110170	EXF-IMRT/SITTS	Exafix IMRT Baseplate, fit for Siemens TTS couchtop
110180	EXF-IMRT/EIMI	Exafix IMRT Baseplate, fit for Elekta MI-beam evo couchtop
110190	EXF-IMRT/ELPR	Exafix IMRT Baseplate, fit for Elekta Precise/C-arm couchtop
110280	EXF-IMRT/CF	Exafix IMRT Baseplate, fit on flat table top, carbon fiber, 2p
110290	EXF-IMRT/GF	Exafix IMRT Baseplate, fit on flat table top, glass fiber, 3p

Contact MACROMEDICS or your local distributor for detailed information or an on site product demo.



macromedics

SOLUTIONS IN RADIOTHERAPY

MACROMEDICS BV
Kouwe Hoek 18
2741 PX Waddinxveen
The Netherlands

Tel: +31(0)182 389777
Fax: +31(0)182 389778
info@macromedics.com
www.macromedics.com