

MACROMEDICS®

Head Supports™

The MACROMEDICS® 3-D shaped head supports provide stable and comfortable support during treatment in the head- and neck region. The 3-D shaped head supports are made of a low density foam and are covered with an easy to clean coating. The head supports are colour coded and come in a set of 3. Each head support will provide a specific, preferred position to create the most optimal approach for treatment.

A comfortable wide head support is available for all MACROMEDICS head- and neck, breast, lung and thorax treatment devices. Also available is a special set of pediatric head supports.

For situations where extra stability and reproducibility is desired MACROMEDICS developed a special low attenuation set of 3-D shaped carbon fiber head supports. The head supports consists of a very thin layer of carbon fiber combined with a low density foam part with precise fit in all MACROMEDICS head- and neck baseplates. For special treatment requirements low attenuation carbon fiber Elevation Blocks and Wedges in various heights and angles are available.

Contact MACROMEDICS or your local distributor for detailed information or an on site product demo.


macromedics

SOLUTIONS IN RADIOTHERAPY

MACROMEDICS BV
 Kouwe Hoek 18
 2741 PX Waddinxveen
 The Netherlands

Tel: +31(0)182 389777
 Fax: +31(0)182 389778
 info@macromedics.com
 www.macromedics.com