



# NIKE PLASTIC RADIATION GLASSES

## MODEL RG-NI-BRZN



RG-NI-BRZN-EV0758-049  
*(Black Blue)*



RG-NI-BRZN-EV0571-001  
*(Black)*

### DESCRIPTION

The Nike Brazen frame combines streamlined style, comfort and lightweight durability. It features a nylon frame, metal accents and a ventilated nose bridge for reduced fogging. This product offers our standard high quality, distortion-free SF-6 Schott Glass Radiation-reducing lenses with .75mm lead equivalency. The Model RG-NI-BRZN fits best for Large head sizes.

### FEATURES

- *0.75mm Lead Equivalency Schott SF6 Radiation Protective Lenses*
- *Adjustable Ventilated Nose*
- *Secure Wrap Temple Arms*

### SPECS

*Frame Size: 60-16-135*  
*Weight: 80g*



Phillips Safety Products  
123 Lincoln Blvd.  
Middlesex, NJ 08846

Contact Us  
1-866-575-1307  
info@phillips-safety.com