

# Protective floor screen

Ref. PM0 - PM1 - PM2 - PM3 - PM4 - PM5



PM1



PM2



PM3



PM4



PM5

## Principle

The Lemer Pax mobile screens providing 2 mm lead protection are realised in flame resistant materials with contemporary design.

The 2 mm lead equivalent lead glass window assure an homogenous protection with the lead lined opaque panel finished with a white pearl durable surface. The entire screens are mounted on heavy duty steel legs and pivoting ball bearings casters with brakes.

Reference	06100000	06100004	06100014	06100016	06100018	06100019
Dimensions	PM0	PM1	PM2	PM3	PM4	PM5
Front panel	750 * 2000	750 * 2000	800 * 2000	950 * 1170/2000	730 * 2000	750 * 2000
Lead glass	300 * 200	400 * 600	800 * 800	800 * 830	500 * 800	500 * 800
Side protection						250 * 1710