





# One Cradle, Three Solutions

DSPS, SSPS and DSPS Prone



DSPS® Prone



SSPS®



DSPS®

## Dosimetric properties of the Cradle

The performance specifications are for reference only, values may vary based on the set-up and the equipment used.

Product Type	Energy
	6 MV
DSPS CF Cradle Bottom Plate	0.6%
DSPS CF Cradle Rim	0.6 - 0.9%
DSPS CF Cradle vertical columns	1.1 - 2.1%

Weight (approx):	0.2 kg / 0.45 lbs
Primary Materials:	Carbon Fiber