

BreastBoard SX™ and LX™

Breast and Thorax

The MacroMedics® Breast treatment line includes the BreastBoard SX and LX. These boards have a straightforward design for quick, easy and efficient patient set-up. Utilizing a coreless, single carbon fiber layer for low attenuation, the BreastBoards provide an IGRT treatment zone with maximal homogeneity, rigidity and stability.



BreastBoard SX

The BreastBoard SX is the elementary model of the MacroMedics® breast treatment line. It is a cost-effective and modern solution: an ergonomic breast treatment device securing reproducible day to day patient set-up.



BreastBoard LX

The BreastBoard LX is the MacroMedics premium, stand-alone breastboard. The stable, easy to adjust arm supports on the BreastBoard LX address many of the adjustment challenges presented by other breast board arm supports, providing quick, customizable set-ups for patient comfort. The BreastBoard LX further addresses common patient comfort issues by incorporating a contoured lumbar area and an optional bottomstop adapter for shorter patients.



Bottomstop



Bottomstop with
bottomstop adaption

