

Scanflex | *Medical*



Flexi -Keep

a mobile system for hanging and storage puts wheels
on your protective clothing

Flexi -Keep


Store your radiation protective clothing
easily and conveniently




Art.nr 10.005



Art.nr 10.005D



A hanger for skirt
Art.nr 10.110 G



Art.nr 10.010



Art.nr 10.010D



Mod.nr	Model	Number of hangers	Height	Width
10.005	Flexi-Keep, standard	5	1460 mm	530 mm
10.005 D	Flexi-Keep, standard	10	1460 mm	530 mm
10.005 T	Rail	1-5		530 mm
10.010	Flexi-Keep, standard	10	1460 mm	880 mm
10.010 D	Flexi-Keep, standard	20	1460 mm	880 mm
10.110 G	Hanger for skirt	1		

Scanflex | *Medical*