

High sensitivity hand-held metal detector

MRI compatible



0,58 kg

ADVANTAGES

- Detects magnetic and non-magnetic metals

Lemer Pax Extras

- Very long battery life, charger
- Visual, audible and vibrating signal
- Exclusive ergonomic shape

CHARACTERISTICS

- Automatic calibration
- Built-in battery
- Shipping dimensions:
L. 440 x H. 160 x D. 270 mm
- Gross weight with charger: 2.5 kg
- Item Code:
 - 02050008 Normal sensitivity
 - 02050010 High sensitivity
 - 02050013 Normal sensitivity with inductive battery charger
 - 02050014 High sensitivity with inductive battery charger for serial connection
 - 02050009 Inductive battery charger

